

## FAQ'S

### WHAT IS SKIN NEEDLING?

Skin Needling is a non-invasive, clinical treatment that uses microscopic needles to treat a variety of skin conditions ranging from ageing and pigmentation to acne, scarring and stretch marks.

### DOES IT HURT?

skin-inject™ DNC Derma Roller Homecare treatments require no downtime and have no discomfort. A numbing cream will be applied to increase comfort when receiving a skin-inject™ DNC Derma Roller clinical treatment. Patients can resume normal activities almost immediately, but may experience some minor redness for a day or two.

### WHO IS IT FOR?

skin-inject™ DNC Derma Roller has been designed for any person experiencing such skin conditions including ageing; pigmentation; scarring; acne; enlarged pores or stretch marks.

### WHAT RESULTS CAN I SEE?

Results may be seen even after one treatment. Patients report a firmer, even-toned, more revitalised skin with a visible decrease in the appearance of wrinkles, pigmentation, sun damage, acne, enlarged pores and scarring.

### HOW MANY TREATMENTS WILL I NEED?

Depending on the appropriate course prescribed by your skin-inject™ DNC Derma Roller practitioner, you may experience anywhere between 3-8 treatments with 6-8 week intervals between each visit.

### HOW MUCH DOES IT COST?

Your treatment prescriber will provide you with a full recommended treatment course and costing upon consultation.

## RECOVERY

The needle length and gauge of all skin-inject™ DNC Derma Roller devices have been specifically engineered to ensure the shortest and most comfortable recover time possible, whilst still ensuring maximum results and treatment longevity.

There is no downtime when using the skin-inject™ DNC Derma Roller Homecare device.

Immediately following in-clinic treatments, the skin may appear temporarily reddened, similar to that of a mild to light sunburn. Whilst some patients have reported that the skin remains red and slightly swollen for 2-4 days, there will be no scabbing or lesions present.

## CONSIDERATIONS/CONTRAINDICATIONS

As skin-inject™ DNC Derma Roller is a cosmedical device designed to deliver dramatic and long-lasting results, a thorough consultation by a qualified skin-inject™ DNC Derma Roller practitioner will ensure complete safety and compatibility. Certain medications and health conditions may affect skin-inject™ DNC Derma Roller clinical treatments. It is recommended to avoid any active skin care or treatment for a minimum of 4 days post a skin-inject™ DNC Derma Roller clinical appointment.

## COMFORT

As skin-inject™ DNC Derma Roller clinical treatments are performed at a microscopic level within the skin, there is minimal discomfort. All treatments are enhanced with the use of a numbing cream to provide further care and patient satisfaction.

skin-inject™ DNC Derma Roller Homecare devices provide a completely comfortable experience and do not require use of any numbing cream.

## RESULTS

### VISIBLE SKIN REJUVENATION

### SOFTER AND SMOOTHER FEELING SKIN

### REDUCED APPEARANCE OF:

- AGEING
- UNEVEN PIGMENTATION
- SCARRING
- STRETCH MARKS
- ENLARGED PORES
- BREAKOUT
- CAPILLARIES



SKININJECT.COM

INFO@SKININJECT.COM






skin-inject™ DNC Derma Roller utilises 600 gamma-sterilised, surgical-grade stainless steel micro needles with patented Disc Needle Care technology to rejuvenate, regenerate and revitalise the skin's own natural healing response. Ten independently rotating discs with an advanced silicone suspension system provide the most comfortable and beneficial treatment possible, whilst reducing procedure downtime.

skin-inject™ DNC Derma Roller offers 4 different treatment lengths providing results that may be performed personally at home or under the care of a skin-inject™ DNC Derma Roller practitioner. skin-inject™ DNC Derma Roller's exclusive design holds 4 patents and is internationally certified and accredited with FDA, KFDA, CE,ISO and GMP approval.

skin-inject™ DNC Derma Rollers provide the fastest, safest and most comfortable skin needling experience.



## CONCEPT

Skin needling is originally derived from ancient acupuncture techniques, and was developed in the early 1990s by a team of plastic surgeons, medical engineers and dermatologists. skin-inject™ DNC Derma Roller is a non-invasive procedure that utilises small rolling discs containing microscopic needles between 0.25mm to 1.5mm depending on the model used. Tiny puncture channels are created within the epidermis and dermis to encourage an organic renewal and regenerating response. There are no injections, incisions, stitches or implants made or inserted.

## BENEFITS

skin-inject™ DNC Derma Roller clinical treatments unlike Laser resurfacing may be performed on all skin colours and backgrounds with minimum downtime. For the best treatment plan, it is always recommended to consult with a skin-inject™ DNC Derma Roller Practitioner.

**SKIN-INJECT™ DNC DERMA ROLLER DEVICES HAVE BEEN DESIGNED TO ADDRESS THE VISUAL MANAGEMENT OF:**

- SKIN REJUVENATION
- TOPICAL PRODUCT INFUSION
- AGEING / WRINKLES
- PIGMENTATION / UV DAMAGE
- ACNE / PIMPLES / BLACKHEADS
- SCARRING
- ENLARGED PORES / SEBORRHOEA
- STRETCH MARKS

skin-inject™ DNC Derma Roller treatments also assist in the increased appearance of a more refreshed, revived and awakened complexion.

There are 4 different sizes of skin-inject™ DNC Derma Roller devices. A qualified practitioner will prescribe the appropriate course of treatments.

## TREATMENTS

**AT HOME:** skin-inject™ DNC Derma Roller Homecare can be used 2-3 times a week at home. Treatments can be performed in as little as 3 minutes and the device will last up to 6 months before a replacement is required.

**CLINICAL:** After a thorough skin analysis, a skin-inject™ DNC Derma Roller practitioner will prescribe a personal course of needling treatments appropriate to the skin conditions being treated. Ideally treatments are performed every 6-8 weeks to ensure maximum results. Each clinical appointment will last approximately 45-60 minutes.

## SCIENCE

As the surgical grade stainless steel needles are gently rolled over the skin, a natural skin reaction is created where circulation is increased and tiny electrical signals are triggered.

This creates a biological chain of events where essential minerals, proteins, building blocks and collagen are discharged and regenerated. After a course of treatments, the skin is visibly revitalised, strengthened and replenished.

## TREATMENT PROFILE GUIDE

TREATMENT	0.25mm Homecare	0.5mm Clinical	1.0mm Clinical	1.5mm Clinical
Collagen Induction Therapy (CIT)	✓	✓	✓	✓
Product Infusion / Skin Stimulation	✓	✓	✓	✓
Anti-Ageing / Wrinkles	✓	✓	✓	✓
Pigmentation / Sun Damage	✓	✓	✓	✓
Rosacea		✓		
Acne / Pimples / Blackheads	✓	✓	✓	
Scarring			✓	✓
Enlarged Pores / Seborrhoea			✓	✓
Stretch Marks				✓
Capillaries / Angiogenesis		✓	✓	✓
Hair Loss / Thinning Hair	✓			